

Product datasheet for AR50772PU-N

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436

Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OCT4 (1-265, His-tag) Human Protein

Product data:

Product Type: Recombinant Proteins

Description: OCT4 (1-265, His-tag) human recombinant protein, 0.5 mg

Species: Human
Expression Host: E. coli

Expression cDNA Clone MGSSHHHHHH SSGLVPRGSH MHFYRLFLGA TRRFLNPEWK GEIDNWCVYV LTSLLPFKIQ

or AA Sequence: SQDIKALQKE LEQFAKLLKQ KRITLGYTQA DVGLTLGVLF GKVFSQTTIC RFEALQLSFK NMCKLRPLLQ

KWVEEADNNE NLQEICKAET LVQARKRKRT SIENRVRGNL ENLFLQCPKP TLQQISHIAQ QLGLEKDVVR VWFCNRRQKG KRSSSDYAQR EDFEAAGSPF SGGPVSFPLA PGPHFGTPGY

GSPHFTALYS SVPFPEGEAF PPVSVTTLGS PMHSN

Tag: His-tag
Predicted MW: 32.2 kDa
Concentration: lot specific

Purity: >85% by SDS - PAGE

Buffer: Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

Preparation: Liquid purified protein

Protein Description: Recombinant human OCT4 protein, fused to His-tag at N-terminus, was expressed in E.coli.

Storage: Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid

repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeg: NP 001167002

 Locus ID:
 5460

 UniProt ID:
 M1S623

 Cytogenetics:
 6p21.33

Synonyms: Oct-3; Oct-4; OCT3; OCT4; OTF-3; OTF4






Summary:

This gene encodes a transcription factor containing a POU homeodomain that plays a key role in embryonic development and stem cell pluripotency. Aberrant expression of this gene in adult tissues is associated with tumorigenesis. This gene can participate in a translocation with the Ewing's sarcoma gene on chromosome 21, which also leads to tumor formation. Alternative splicing, as well as usage of alternative AUG and non-AUG translation initiation codons, results in multiple isoforms. One of the AUG start codons is polymorphic in human populations. Related pseudogenes have been identified on chromosomes 1, 3, 8, 10, and 12. [provided by RefSeq, Oct 2013]

Protein Families:

Adult stem cells, Cancer stem cells, Embryonic stem cells, Induced pluripotent stem cells, Stem cell - Pluripotency, Transcription Factors

Product images:



15% SDS-PAGE (3ug)