

Product datasheet for **AR50795PU-S**

CDH11 (54-617, His-tag) Human Protein


Product data:

Product Type:	Recombinant Proteins
Description:	CDH11 (54-617, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Tag:	His-tag
Predicted MW:	66.2 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH8.0) containing 10% glycerol 0.1M NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant human CDH11 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Locus ID:	109864281
Synonyms:	RNA5-8N2



[View online »](#)

Product images:



15% SDS-PAGE (3ug)