

Product datasheet for **AR50816PU-S**

NRIP3 (1-241, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	NRIP3 (1-241, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMFYSGLL TEGGRKETDM REAASLRQQR RMKQAVQFIH KDSADLLPLD GLKKGSSKD MQPHNILQRR LMETNLSKLR SGPRVPWASK TNKLNQAKSE GLKKSEEDDM ILVSCQCAGK DVKALVDTGC LYNLISLACV DRLGLKEHVK SHKHEGEKLS LPRHLKVGQ IEHLVITLGS LRLDCPAAW DDNEKNLSLG LQTLRSLKCI INLDKHRLIM GKTDKEEIPF VETVSLNEDN TSEA
Tag:	His-tag
Predicted MW:	29.4 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human NRIP3 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_065696
Locus ID:	56675
UniProt ID:	Q9NQ35
Cytogenetics:	11p15.4
Synonyms:	C11orf14; NY-SAR-105



[View online »](#)

Summary: NRIP3, also known as C11orf14, is a 241 amino acid protein that is encoded by a gene which maps to human chromosome 11.

Protein Families: Druggable Genome

Product images:

