

## Product datasheet for **AR50888PU-S**

### CSDC2 (1-153, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	CSDC2 (1-153, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMTSESTS PPVVPPLHSP KSPVWPTFPF HREGSRVWER GGVPPRDLPs PLPTKRTRTY SATARASAGP VFKGVCKQFS RSQGHGFITP ENGSEDIFVH VSDIEGEYVP VEGDEVYKYM CPIPKNQKF QAVEVLTQL APHTPHETWS GQVVGs
Tag:	His-tag
Predicted MW:	19.2 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol 0.1M NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant human CSDC2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_055275</a>
Locus ID:	27254
UniProt ID:	<a href="#">Q9Y534</a>
Cytogenetics:	22q13.2
Synonyms:	dj347H13.2; PIPPIN



[View online »](#)

**Summary:** RNA-binding factor which binds specifically to the very 3' UTR ends of both histone H1 and H3.3 mRNAs, encompassing the polyadenylation signal. Might play a central role in the negative regulation of histone variant synthesis in the developing brain (By similarity). [UniProtKB/Swiss-Prot Function]

**Protein Families:** Transcription Factors

**Product images:**

