

Product datasheet for **AR50912PU-S**

OCIAD2 (1-154, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	OCIAD2 (1-154, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH</u> <u>SSGLVPRGSH</u> <u>MGSMASASAR</u> GNQDKDAHFP PPSKQSLFC PKSKLHIHRA EISKIMRECQ EESFWKRALP FSLVSMLVTQ GLVYQGYLAA NSRFGSLPKV ALAGLLGFGL GKVSYIGVCQ SKFHFFEDQL RGAGFGPQHN RHCLLTCEEC KIKHGLSEKG DSQPSAS
Tag:	His-tag
Predicted MW:	19.3 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.4M UREA, 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant Human OCIAD2 protein, fused to His-tag at N-terminus, was expressed in E.coli.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_001014446</u>
Locus ID:	132299
UniProt ID:	<u>Q56VL3</u>
Cytogenetics:	4p11
Summary:	OCIAD2 is a protein that localizes to endosomes and contains one OCIA (Ovarian carcinoma immunoreactive antigen) domain.



[View online »](#)

Product images:

