

## Product datasheet for **AR50933PU-S**

### **RABL5 (1-185, His-tag) Human Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	RABL5 (1-185, His-tag) human protein, 20 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMLKAKIL FVGPCESGKT VLANFLTESS DITEYSPTQG VRILEFENPH VTSNNKGTGC EFELWDCGGD AKFESCWPAL MKDAHGVVIV FNADIPSHRK EMEMWYSCFV QQPSLQDTQC MLIHHKPGS GDDKGSLSLS PPLNKLKLVH SNLEDDPEEI RMEFIKYLKS IINSMSESRD REEMSIMT
Tag:	His-tag
Predicted MW:	23.2 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_001124292</a>
Locus ID:	64792
UniProt ID:	<a href="#">Q9H7X7</a>
Cytogenetics:	7q22.1
Synonyms:	FAP9; RABL5
Summary:	Small GTPase-like component of the intraflagellar transport (IFT) complex B. [UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome



[View online »](#)

## Product images:

