

Product datasheet for **AR50944PU-S**

CCDC90B (43-230, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	CCDC90B (43-230, His-tag) human protein, 20 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSYDRRPV DITPLEQRKL TFDTHALVQD LETHGFDKTQ AETIVSALTA LSNVSLDTIY KEMVTQAQQE ITVQQLMAHL DAIRKDMVIL EKSEFANLRA ENEKMKIELD QVKQQLMHET SRIRADNKLD INLERSRVTD MFTDQEKQLM ETTTEFTKKD TQTKSISET SNKIDAEIAS LKTLMESNKL E
Tag:	His-tag
Predicted MW:	24.1 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001273045
Locus ID:	60492
UniProt ID:	Q9GZT6
Cytogenetics:	11q14.1
Synonyms:	MDS011; MDS025
Summary:	CCDC90B belongs to the CCDC90 family. This protein is characterized by the presence of a domain of unknown function, DUF1640. This domain is a characteristic of the entire protein with the exception of the first twenty-three amino acid residues which are



[View online »](#)

Protein Families: Transmembrane

Product images:

