

## Product datasheet for **AR50946PU-S**

### ANKRD54 / LIAR (1-300, His-tag) Human Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	ANKRD54 / LIAR (1-300, His-tag) human recombinant protein, 20 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH</u> <u>SSGLVPRGSH</u> MGSMAAAAGD ADDEPRSGHS SSEGECAVAP EPLTDAEGLF SFADFGSALG GGGAGLSGRA SGGAQSPRY LHVLWQQDAE PRDELCKIP AGRLLRAARP HRRLGPTGKE VHALKRLRDS ANANDVETVQ QLLEDGADPC AADDKGR TAL HFASCNGNDQ IVQLLLDHGA DPNQRDGLGN TPLHLA ACTN HVPVITLLR GGARVDALDR AGRTPLHLAK SKLNILQEGH AQCLEAVRLE VKQIIHMLRE YLERLGQHEQ RERLDDLCTR LQMTSTKEQV DEVTDLLASF TSLSLQMQSM EKR
Tag:	His-tag
Predicted MW:	34.9 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human ANKRD54 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_001336782</u>
Locus ID:	129138
Cytogenetics:	22q13.1
Synonyms:	LIAR



[View online »](#)

**Summary:** Plays an important role in regulating intracellular signaling events associated with erythroid terminal differentiation.[UniProtKB/Swiss-Prot Function]

**Product images:**



15% SDS-PAGE (3ug)