

## Product datasheet for **AR50947PU-N**

### COMMD7 (1-200, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	COMMD7 (1-200, His-tag) human protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMGRLHCT EDPVPEAVGG DMQQLNQLGA QQFSALTEVL FHFLTEPKVEV ERFLAQLSEF ATTNQISLGS LRSIVKSLLL VPNGALKKSL TAKQVQADFI TLGLSEEKAT YFSEKWKQNA PTLARWAIGQ TLMINQLIDM EWKFGVTS GS SELEKVGSI F LQLKLVVKKG NQTENVYIEL TLPQFYSFLH EMERVTSME CFC
Tag:	His-tag
Predicted MW:	24.9 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_001092809</a>
Locus ID:	149951
UniProt ID:	<a href="#">Q86VX2</a>
Cytogenetics:	20q11.21
Synonyms:	C20orf92; dj1085F17.3
Summary:	May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes (PubMed:21778237). Associates with the NF-kappa-B complex and suppresses its transcriptional activity (PubMed:15799966).[UniProtKB/Swiss-Prot Function]



[View online »](#)

## Product images:

