

Product datasheet for AR50959PU-S

BUD31 / EDG2 (1-144, His-tag) Human Protein

Product data:

OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Recombinant Proteins
Description:	BUD31 / EDG2 (1-144, His-tag) human protein, 20 μg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMPKVKRS RKAPPDGWEL IEPTLDELDQ KMREAETEPH EGKRKVESLW PIFRIHHQKT RYIFDLFYKR KAISRELYEY CIKEGYADKN LIAKWKKQGY ENLCCLRCIQ TRDTNFGTNC ICRVPKSKLE VGRIIECTHC GCRGCSG
Tag:	His-tag
Predicted MW:	19.4 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% glycerol, 2 mM DTT
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP 003901</u>
Locus ID:	8896
UniProt ID:	<u>P41223</u>
Cytogenetics:	7q22.1
Synonyms:	Cwc14; EDG-2; EDG2; fSAP17; G10; YCR063W
Summary:	Involved in the pre-mRNA splicing process (PubMed:28502770, PubMed:28076346). May play a role as regulator of AR transcriptional activity; may increase AR transcriptional activity (PubMed:25091737).[UniProtKB/Swiss-Prot Function]
Protein Families:	Transcription Factors



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US Protein Pathways: Spliceosome

Product images:



15% SDS-PAGE (3ug)

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US