

Product datasheet for **AR50967PU-S**

Stathmin-4 / STMN4 (1-216, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Stathmin-4 / STMN4 (1-216, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMTLAAYK EKMKELPLVS LFCSCFLADP LNKSSYKYEG WCGRQCRRKD ESQRKDSADW RERRAQADTV DLNWCVISDM EVIELNKCTS GQSFEVILKP PSFDGVPEFN ASLPRRRDPS LEEIQKKLEA AEERRKYQEA ELLKHLAEKR EHEREVIQKA IEENNNFIKM AKEKLAQKME SNKENREAHL AAMLRLQEK DKHAEVRKN KELKEEASR
Tag:	His-tag
Predicted MW:	27.8 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 20% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human STMN4 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001269982
Locus ID:	81551
UniProt ID:	Q9H169
Cytogenetics:	8p21.2
Synonyms:	RB3
Summary:	Exhibits microtubule-destabilizing activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

