

## Product datasheet for **AR50984PU-N**

### ECRG-2 (20-85, His-tag) Human Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	ECRG-2 (20-85, His-tag) human protein, 0.25 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSSEAASLS PKKVDCSIYK KYPVVAIPCP ITYLPVCGSD YITYGNECHL CTESLKSNGR VQFLHDGSC
Tag:	His-tag
Predicted MW:	9.6 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% glycerol, 2 mM DTT
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_115955</a>
Locus ID:	84651
UniProt ID:	<a href="#">P58062</a>
Cytogenetics:	5q32
Synonyms:	ECG2; ECRG2
Summary:	Probable serine protease inhibitor.[UniProtKB/Swiss-Prot Function]
Protein Families:	Secreted Protein



[View online »](#)

Product images:



15% SDS-PAGE (3ug)