

## **Product datasheet for AR51033PU-S**

## 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

## RPS24 (1-130, His-tag) Human Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** RPS24 (1-130, His-tag) human protein, 20 μg

Species: Human
Expression Host: E. coli

**Expression cDNA Clone** MGSSHHHHHH SSGLVPRGSH MGSMNDTVTI RTRKFMTNRL LQRKQMVIDV LHPGKATVPK

or AA Sequence: TEIREKLAKM YKTTPDVIFV FGFRTHFGGG KTTGFGMIYD SLDYAKKNEP KHRLARHGLY

EKKKTSRKQR KERKNRMKKV RGTAKANVGA GKK

Tag: His-tag
Predicted MW: 17.5 kDa
Concentration: lot specific

Purity: >85% by SDS - PAGE

**Buffer:** Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% glycerol, 2 mM DTT

**Preparation:** Liquid purified protein

Storage: Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid

repeated freezing and thawing.

**Stability:** Shelf life: one year from despatch.

**RefSeq:** NP 001017

 Locus ID:
 6229

 UniProt ID:
 P62847

 Cytogenetics:
 10q22.3

Synonyms: DBA3; eS24; S24





**Summary:** 

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S24E family of ribosomal proteins. It is located in the cytoplasm. Multiple transcript variants encoding different isoforms have been found for this gene. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Mutations in this gene result in Diamond-Blackfan anemia. [provided by RefSeq, Nov 2008]

Protein Pathways:

Ribosome

## **Product images:**

