

Product datasheet for **AR51070PU-N**

Layilin (22-235, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Layilin (22-235, His-tag) human protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSATGRLLS ASDLDLRGGQ PVCRRGGTQRP CYKVIYFHDT SRRLNFEEAK EACRRDGGQL VSIESEDEQK LIEKFENLL PSDGDFWIGL RRREEKQSNSTACQDLYAWT DGSISQFRNW YVDEPSCGSE VCVVMYHQPS APAGIGGPYM FQWNDDRCNM KNNFICKYSD EKPAVPSREA EGEETELTTP VLPEETQEED AKKTFKESRE AA
Tag:	His-tag
Predicted MW:	26.7 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001245319
Locus ID:	143903
UniProt ID:	Q6UX15
Cytogenetics:	11q23.1
Summary:	Receptor for hyaluronate.[UniProtKB/Swiss-Prot Function]
Protein Families:	Transmembrane



[View online »](#)

Product images:

