

Product datasheet for **AR51072PU-S**

ATP synthase subunit d (1-161, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	ATP synthase subunit d (1-161, His-tag) human protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMAGRKLA LKTIDWVFA EIIPQNQKAI ASSLKSWNET LTSRLAALPE NPPAIDWAYY KANVAKAGLV DDFEKKFNAL KVPVPEDKYT AQVDAEEKED VKSCAEWVSL SKARIVEYEK EMEKMKNLIP FDQMTIEDLN EAFPETKLDK KKYPYWPHQP IENL
Tag:	His-tag
Predicted MW:	20.9 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% glycerol
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001003785
Locus ID:	10476
UniProt ID:	O75947
Cytogenetics:	17q25.1
Synonyms:	APT5H; ATP5H; ATPQ



[View online »](#)

Summary:

Mitochondrial ATP synthase catalyzes ATP synthesis, utilizing an electrochemical gradient of protons across the inner membrane during oxidative phosphorylation. It is composed of two linked multi-subunit complexes: the soluble catalytic core, F1, and the membrane-spanning component, Fo, which comprises the proton channel. The F1 complex consists of 5 different subunits (alpha, beta, gamma, delta, and epsilon) assembled in a ratio of 3 alpha, 3 beta, and a single representative of the other 3. The Fo seems to have nine subunits (a, b, c, d, e, f, g, F6 and 8). This gene encodes the d subunit of the Fo complex. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene. In addition, three pseudogenes are located on chromosomes 9, 12 and 15. [provided by RefSeq, Jun 2010]

Protein Pathways:

Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease

Product images: