

Product datasheet for **AR51078PU-S**

CEACAM21 (35-240, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	CEACAM21 (35-240, His-tag) human protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSWLFIASA PFEVAEGENV HLSVVYLPEN LYSYGWYK GK TVEPNQLIAA YVIDTHVRTP GPAYSGRETI SPSGDLHFQN VTLED TGYYN LQV TYRNSQI EQASHHLRVY ESVAQPSIQA SSTTVTEKGS VVLTCHTNNT GTSFQWIFNN QRLQVTKRMK LSWFNHVLT I DPIRQEDAGE YQCEVSNPVS SNRSDPLKLT VKSDDNTLG
Tag:	His-tag
Predicted MW:	25 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% glycerol
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001091976
Locus ID:	90273
UniProt ID:	Q3KPI0 , A0A0G2JSC8
Cytogenetics:	19q13.2
Synonyms:	CEACAM3; R29124_1
Summary:	CEACAM21 is a member of the family of carcinoembryonic antigen-related cell adhesion molecules (CEACAMs), which are used by several bacterial pathogens to bind and invade host cells. The transmembrane protein directs phagocytosis of several bacterial spec



[View online »](#)

Protein Families: Transmembrane

Product images:

