

Product datasheet for **AR51139PU-S**

FAM107B (1-131, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	FAM107B (1-131, His-tag) human protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMAEPDYI EDDNPELIRP QKLINPVKTS RNHQDLHREL LMNQKRGLAP QNKPELQKVM EKRKRQDQVIK QKEEEAQKKK SDLEIELLKR QQKLEQLELE KQKLQEEQEN APEFVKVKG N LRRTGQEVAQ AQES
Tag:	His-tag
Predicted MW:	17.9 kDa
Concentration:	lot specific
Purity:	>90 % by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001269624
Locus ID:	83641
UniProt ID:	Q9H098
Cytogenetics:	10p13
Synonyms:	C10orf45; HITS
Summary:	FAM107B protein, also known as FAM107B, is expressed in most tissues in the human body without there being a high frequency in any one tissue. It is found in all stages of human development. Recombinant human FAM107B protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques



[View online »](#)

Product images:

