

Product datasheet for **AR51140PU-S**

NRBF2 (45-287, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	NRBF2 (45-287, His-tag) human recombinant protein, 20 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMKLTQSE QAHLSLELQR DSHMKQLLI QERWKRAQRE ERLKAQQNTD KDAAAHLQTS HKPSAEDAEG QSPLSQKYSP STEKCLPEIQ GIFDRDPDTL LYLLQQKSEP AEPCIGSKAP KDDKTIIEEQ ATKIADLKRH VEFLLVAENER LRKENKQLKA EKARLLKGPI EKELDVDADF VETSELWSLP PHAETATASS TWQKFAANTG KAKDIPINL PPLDFPSPPEL PLMELSEDIL KGFMMN
Tag:	His-tag
Predicted MW:	29.9 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human NRBF2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001269334
Locus ID:	29982
UniProt ID:	Q96F24
Cytogenetics:	10q21.3
Synonyms:	COPR; COPR1; COPR2; NRBF-2



[View online »](#)

Summary: May modulate transcriptional activation by target nuclear receptors. Can act as transcriptional activator (in vitro).[UniProtKB/Swiss-Prot Function]

Protein Families: Druggable Genome

Product images:

