

## Product datasheet for **AR51262PU-S**

### **FAM84B (1-310, His-tag) Human Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	FAM84B (1-310, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSSLVPRGSH MGSMGNQVEK LTHLSYKEVP TADPTGVDRD DGPRIGVSYI FSNDDDEVDP QPPPQGPDDG GLPDGGDGGP PPQPQPYDPR LHEVECSVFY RDECIYQKSF APGSAALSTY TPENLLNKCK PGDLVEFVSQ AQYPHWAVYV GNFQVHLHR LEVINSFLTD ASQGRRRGRVV NDLYRYKPLS SSAVVRNALA HVGAKERELS WRNSESFAAW CRYGKREFKI GGELRIGKQP YRLQIQLSAQ RSHTLEFQSL EDLIMEKRRN DQIGRAAVLQ ELATHLHPAE PEEGDSNVAR TTPPPGRPPA PSSEEDGGEA VAH
Tag:	His-tag
Predicted MW:	36.9 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol 0.15M NaCl, 1 mM DTT, 0.1 mM PMSF
Preparation:	Liquid purified protein
Protein Description:	Recombinant human FAM84B protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_777571</a>
Locus ID:	157638
UniProt ID:	<a href="#">Q96KN1</a>
Cytogenetics:	8q24.21



[View online »](#)

**Synonyms:** BCMP101; FAM84B; NSE2

**Summary:** FAM84B, as known as Family with sequence similarity 84, member B, coprecipitated with a downstream effector of RAS, CRAF. Binding of FAM83B with CRAF disrupted CRAF/14-3-3 interactions and increased CRAF membrane localization, resulting in elevated MAPK and mammalian target of rapamycin (mTOR) signaling. It is an oncogene and potentially represents a new target for therapeutic intervention.

**Product images:**

