

Product datasheet for **AR51274PU-N**

SEC61B (1-70, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	SEC61B (1-70, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMPGPTPS GTNVGSSGRS PSKAVAAARAA GSTVRQRKNA SCGTRSAGRT TSAGTGGMW R FYTEDSPGLK VGP
Tag:	His-tag
Predicted MW:	9.4 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% glycerol, 2 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human SEC61B protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_006799
Locus ID:	10952
UniProt ID:	P60468
Cytogenetics:	9q22.33



[View online »](#)

Summary:

The Sec61 complex is the central component of the protein translocation apparatus of the endoplasmic reticulum (ER) membrane. Oligomers of the Sec61 complex form a transmembrane channel where proteins are translocated across and integrated into the ER membrane. This complex consists of three membrane proteins- alpha, beta, and gamma. This gene encodes the beta-subunit protein. The Sec61 subunits are also observed in the post-ER compartment, suggesting that these proteins can escape the ER and recycle back. There is evidence for multiple polyadenylated sites for this transcript. [provided by RefSeq, Jul 2008]

Protein Families:

Transmembrane

Protein Pathways:

Vibrio cholerae infection

Product images:
