

## Product datasheet for **AR51326PU-S**

### LZIC (1-190, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	LZIC (1-190, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMASRGKT ETSKCLKQNLE EQLDRMQQL QDLEECREEL DTDEYEETKK ETLEQLSEFN DSLKKIMSGN MTLVDELSGM QLAIQAAISQ AFKTPEVIRL FAKKQPGQLR TRLAEMDRDL MVGKLERDLY TQQKVEILTA LRKLGEKLT A DDEAFLSANA GAILSQFEKV STDLGSGDKI LALASFEVEK TKK
Tag:	His-tag
Predicted MW:	23.9 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human LZIC protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_001303902</a>
Locus ID:	84328
UniProt ID:	<a href="#">Q8WZA0</a> , <a href="#">A0A024R4I7</a>
Cytogenetics:	1p36.22
Summary:	LZIC belongs to the CTNNBIP1 family. This protein is ubiquitously expressed, with highest levels in kidney and up-regulated in several cases of gastric cancers. It does not interact with CTNNB1.



[View online »](#)

## Product images:

