

## Product datasheet for **AR51354PU-S**

### EAF1 (182-268, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	EAF1 (182-268, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSDIIEQMS SSSGSSSSDS ESSSGSDDDS SSSGGEDNGP ASPPQPSHQQ PYNSRPAVAN GTSRPQGSNQ LMNTRLRNDLQ LSESGSDSDD
Tag:	His-tag
Predicted MW:	11.3 kDa
Concentration:	lot specific
Purity:	>90 % by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human EAF1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_149074</a>
Locus ID:	85403
UniProt ID:	<a href="#">Q96JC9</a>
Cytogenetics:	3p25.1
Summary:	Acts as a transcriptional transactivator of ELL and ELL2 elongation activities.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

## Product images:

