

Product datasheet for **AR51361PU-S**

Complement factor B (26-259, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Complement factor B (His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSQ SHMTPWSLAR PQGSCSLEGV EIKGGSFRLL QEQQALEYVC PSGFYPYPVQ TRTCRSTGSW STLKTQDQKT VRKAECRAIH CPRPHDFENG EYWPRSPYYN VSDEISFHCY DGYTLRGSAN RTCQVNGRWS GQTAICDNGA GYCSNPGIPI GTRKVGGSQYR LEDSVTYHCS RGLTLRGSQR RTCQEGGSWS GTEPSCQDSF MYDTPQEVAE AFLSSLTETI EGVDAEDGHG PGEQQKR
Tag:	His-tag
Predicted MW:	28.4 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: This purified protein is available in a denatured form, making it less suitable for functional studies. Denatured proteins are better suited for applications like Western Blot (WB) or imaging assays. State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.4M Urea
Preparation:	Liquid purified protein
Protein Description:	Recombinant human CFB protein, fused to His-tag at N-terminus, was expressed in E.coli.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001701
Locus ID:	629
UniProt ID:	P00751
Cytogenetics:	6p21.33
Synonyms:	Properdin factor B, C3/C5 convertase, PBF2, CFB, BF, BFD



[View online »](#)

Summary:

This gene encodes complement factor B, a component of the alternative pathway of complement activation. Factor B circulates in the blood as a single chain polypeptide. Upon activation of the alternative pathway, it is cleaved by complement factor D yielding the noncatalytic chain Ba and the catalytic subunit Bb. The active subunit Bb is a serine protease which associates with C3b to form the alternative pathway C3 convertase. Bb is involved in the proliferation of preactivated B lymphocytes, while Ba inhibits their proliferation. This gene localizes to the major histocompatibility complex (MHC) class III region on chromosome 6. This cluster includes several genes involved in regulation of the immune reaction. Polymorphisms in this gene are associated with a reduced risk of age-related macular degeneration. The polyadenylation site of this gene is 421 bp from the 5' end of the gene for complement component 2. [provided by RefSeq, Jul 2008]

Protein Families:

Druggable Genome, Protease, Secreted Protein

Protein Pathways:

Complement and coagulation cascades

Product images: