

# Product datasheet for AR51362PU-N

## MMP-1 (100-469, His-tag) Human Protein

### **Product data:**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Recombinant Proteins
Description:	MMP-1 (100-469, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSFVLTEGN PRWEQTHLTY RIENYTPDLP RADVDHAIEK AFQLWSNVTP LTFTKVSEGQ ADIMISFVRG DHRDNSPFDG PGGNLAHAFQ PGPGIGGDAH FDEDERWTNN FREYNLHRVA AHELGHSLGL SHSTDIGALM YPSYTFSGDV QLAQDDIDGI QAIYGRSQNP VQPIGPQTPK ACDSKLTFDA ITTIRGEVMF FKDRFYMRTN PFYPEVELNF ISVFWPQLPN GLEAAYEFAD RDEVRFFKGN KYWAVQGQNV LHGYPKDIYS SFGFPRTVKH IDAALSEENT GKTYFFVANK YWRYDEYKRS MDPGYPKMIA HDFPGIGHKV DAVFMKDGFF YFFHGTRQYK FDPKTKRILT LQKANSWFNC RKN
Tag:	His-tag
Predicted MW:	45.0 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.4M Urea
Preparation:	Liquid purified protein
Protein Description:	Recombinant human MMP1 protein, fused to His-tag at N-terminus, was expressed in E.coli.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP 001139410</u>
Locus ID:	4312
UniProt ID:	<u>B4DN15</u>
Cytogenetics:	11q22.2
Synonyms:	CLG; CLGN




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	MMP-1 (100-469, His-tag) Human Protein – AR51362PU-N
Summary:	This gene encodes a member of the peptidase M10 family of matrix metalloproteinases (MMPs). Proteins in this family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. The encoded preproprotein is proteolytically processed to generate the mature protease. This secreted protease breaks down the interstitial collagens, including types I, II, and III. The gene is part of a cluster of MMP genes on chromosome 11. Mutations in this gene are associated with chronic obstructive pulmonary disease (COPD). Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed. [provided by RefSeq, Jan 2016]
Protein Familie	es: Druggable Genome, Protease, Secreted Protein

**Protein Pathways:** Bladder cancer, Pathways in cancer, PPAR signaling pathway

# **Product images:**



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US