

Product datasheet for **AR51382PU-N**

BTF3L4 (1-78, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	BTF3L4 (1-78, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMNQEKLA KLQAQVRIGG KGTARRKKKV VHRTATADDK KLQSSLKLA VNNIAGIEES WTVKHQNQKT LMRKMMMFQI L
Tag:	His-tag
Predicted MW:	11.3 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.4M Urea
Preparation:	Liquid purified protein
Protein Description:	Recombinant human BTF3L4 protein isoform 3, fused to His-tag at N-terminus, was expressed in E.coli .
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001129969
Locus ID:	91408
UniProt ID:	Q96K17
Cytogenetics:	1p32.3
Summary:	BTF3L4 belongs to the NAC-beta family and contains 1 NAC-A/B (NAC-alpha/beta) domain.



[View online »](#)

Product images:

