

Product datasheet for **AR51430PU-N**

CSNK1A1 (1-337, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	CSNK1A1 (1-337, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSSLVPRGSH MASSSGSKAE FIVGGKYKLV RKIGSGSFGD IYLAINITNG EEVAVKLESQ KARHPQLLYE SKLYKILQGG VGIPHIRWYG QEKDYNVLVM DLLGPSLEDL FNFCRRFTM KTVLMLADQM ISRIEYVHTK NFIHRDIKPD NFLMGIGRHC NKLFLIDFGL AKKYRDNRTR QHIPYREDKN LTGTARYASI NAHLGIEQSR RDDMESLGYV LMYFNRTSLP WQGLKAATKK QKYEKISEKK MSTPVEVLCK GFPAEFAMYL NYCRGLRFEE APDYMRLRQL FRILFRTLNH QDYTFDWTM LKQKAAQQA SSSGQGQQAQ TPTGKQTDKT KSNMKGF
Tag:	His-tag
Predicted MW:	41 kDa
Concentration:	lot specific
Purity:	>80% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human CSNK1a1 protein, fused to His-tag at N-terminus, was expressed in E.coli.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001020276
Locus ID:	1452
UniProt ID:	P48729 , Q6PJ06 , P48729-2
Cytogenetics:	5q32
Synonyms:	CK1; CK1a; CKIa; HEL-S-77p; HLCDGP1; PRO2975



[View online »](#)

Summary:

Casein kinases are operationally defined by their preferential utilization of acidic proteins such as caseins as substrates. It can phosphorylate a large number of proteins. Participates in Wnt signaling. Phosphorylates CTNNB1 at 'Ser-45'. May phosphorylate PER1 and PER2. May play a role in segregating chromosomes during mitosis (PubMed:11955436, PubMed:1409656, PubMed:18305108). May play a role in keratin cytoskeleton disassembly and thereby, it may regulate epithelial cell migration (PubMed:23902688).[UniProtKB/Swiss-Prot Function]

Protein Families:

Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase

Protein Pathways:

Hedgehog signaling pathway, Wnt signaling pathway

Product images: