

Product datasheet for AR51458PU-S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CD201 / EPCR (18-210, His-tag) Human Protein

Product data:

Product Type: Recombinant Proteins

Description: CD201 / EPCR (18-210, His-tag) human recombinant protein, 0.1 mg

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

MGSSHHHHHH SSGLVPRGSH MGSSQDASDG LQRLHMLQIS YFRDPYHVWY QGNASLGGHL THVLEGPDTN TTIIQLQPLQ EPESWARTQS GLQSYLLQFH GLVRLVHQER TLAFPLTIRC

FLGCELPPEG SRAHVFFEVA VNGSSFVSFR PERALWQADT QVTSGVVTFT LQQLNAYNRT

FLUCELPPEU SKANVFFEVA VNUSSFVSFK PEKALVVQADT QVTSUVVTFT LQQLIV

RYELREFLED TCVQYVQKHI SAENTKGSQT SRSYTS

Tag: His-tag
Predicted MW: 24.4 kDa
Concentration: lot specific

Purity: >80% by SDS - PAGE

Buffer: Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol 0.4M Urea

Preparation: Liquid purified protein

Protein Description: Recombinant human PROCR protein, fused to His-tag at N-terminus, was expressed in E.coli.

Storage: Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid

repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: NP 006395

 Locus ID:
 10544

 UniProt ID:
 Q9UNN8

 Cytogenetics:
 20q11.22

Synonyms: CCCA; CCD41; EPCR






Summary:

The protein encoded by this gene is a receptor for activated protein C, a serine protease activated by and involved in the blood coagulation pathway. The encoded protein is an N-glycosylated type I membrane protein that enhances the activation of protein C. Mutations in this gene have been associated with venous thromboembolism and myocardial infarction, as well as with late fetal loss during pregnancy. The encoded protein may also play a role in malarial infection and has been associated with cancer. [provided by RefSeq, Jul 2013]

Protein Families:

Druggable Genome, ES Cell Differentiation/IPS, Transmembrane

Product images:



15% SDS-PAGE (3ug)