

Product datasheet for **AR51518PU-S**

CNPY4 (22-248, His-tag) Human Protein

Product data:


Product Type:	Recombinant Proteins
Description:	CNPY4 (22-248, His-tag) human recombinant protein, 20 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSEFGMLKE EDDDTERLPS KCEVCKLLST ELQAELSRTG RSREVLGLGQ VLDTGKRKRH VPYSVSETRL EEALENLCEER ILDYSVHAER KGSLRYAKGQ SQTMATLKGL VQKGVKVDLG IPLELWDEPS VEVTYLKKQC ETMLEEFEDI VGDWYFHHQE QPLQNFLCEG HVLPAETAC LQETWTGKEI TDGEEKTEGE EEQEEEEEEE EEEGGDKMTK TGSHPKLDRE DL
Tag:	His-tag
Predicted MW:	28.6 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% glycerol, 5 mM DTT, 2 mM EDTA
Preparation:	Liquid purified protein
Protein Description:	Recombinant human CNPY4, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_689968
Locus ID:	245812
UniProt ID:	Q8N129
Cytogenetics:	7q22.1
Synonyms:	PRAT4B



[View online »](#)

Summary: Plays a role in the regulation of the cell surface expression of TLR4.[UniProtKB/Swiss-Prot Function]

Product images:



15% SDS-PAGE (3ug)