

Product datasheet for **AR51520PU-N**

NAIF1 / C9orf90 (1-327, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	NAIF1 / C9orf90 (1-327, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMVPAKK RKMNFSEREV EIIVEELELK KHLLVNHFNAGVPLAAKSAA WHGILRRVNA VATCRRELPE VKKKWSDLKT EVRRKVAQVR AAVEGGEAPG PTEEDGAGGP GTGGGSGGGG PAVAPVLLTP MQQRICNLLG EATIISLPST TEIHPVALGP SATAAAATVT LTQIPTETTY HTLEEGVVEY CTAEAPPPLP PETPVDMMMAQ HADTSVKPQA LKSRIALNSA KLIQEQRVTN LHVKEIAQHL EQQNDLLQMI RRSQEVQACA QERQAQAMEG TQAALSVDLIQ VLRPMIKDFR RYLQSNTANP APASDPGQVA QNGQPDSIIQ
Tag:	His-tag
Predicted MW:	37.6 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 30% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human NAIF1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_931045
Locus ID:	203245
UniProt ID:	Q69YI7 , A0A024R887
Cytogenetics:	9q34.11




[View online »](#)

Synonyms: bA379C10.2; C9orf90

Summary: Induces apoptosis.[UniProtKB/Swiss-Prot Function]

Product images:



15% SDS-PAGE (3ug)