

Product datasheet for AR51522PU-N

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

TARBP2 (1-345,) Human Protein

Product data:

Product Type: Recombinant Proteins

Description: TARBP2 (1-345,) human recombinant protein, 0.5 mg

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

MLAANPGKTP ISLLQEYGTR IGKTPVYDLL KAEGQAHQPN FTFRVTVGDT SCTGQGPSKK AAKHKAAEVA LKHLKGGSML EPALEDSSSF SPLDSSLPED IPVFTAAAAA TPVPSVVLTR

SPPMELQPPV SPQQSECNPV GALQELVVQK GWRLPEYTVT QESGPAHRKE FTMTCRVERF IEIGSGTSKK LAKRNAAAKM LLRVHTVPLD ARDGNEVEPD DDHFSIGVGS RLDGLRNRGP

GCTWDSLRNS VGEKILSLRS CSLGSLGALG PACCRVLSEL SEEQAFHVSY LDIEELSLSG LCQCLVELST

QPATVCHGSA TTREAARGEA ARRALQYLKI MAGSK

Predicted MW: 36.9 kDa

Concentration: lot specific

Purity: >85% by SDS - PAGE

Buffer: Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol 0.4M Urea

Preparation: Liquid purified protein

Protein Description: Recombinant human TARBP2 protein was expressed in E.coli.

Storage: Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid

repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: NP 004169

Locus ID: 6895

 UniProt ID:
 Q15633

 Cytogenetics:
 12q13.13

Synonyms: LOQS; TRBP; TRBP1; TRBP2






Summary:

HIV-1, the causative agent of acquired immunodeficiency syndrome (AIDS), contains an RNA genome that produces a chromosomally integrated DNA during the replicative cycle. Activation of HIV-1 gene expression by the transactivator Tat is dependent on an RNA regulatory element (TAR) located downstream of the transcription initiation site. The protein encoded by this gene binds between the bulge and the loop of the HIV-1 TAR RNA regulatory element and activates HIV-1 gene expression in synergy with the viral Tat protein. Alternative splicing results in multiple transcript variants encoding different isoforms. This gene also has a pseudogene. [provided by RefSeq, Jul 2008]

Protein Families:

Stem cell - Pluripotency

Product images:



15% SDS-PAGE (3ug)