

Product datasheet for **AR51536PU-S**

COA4 (1-87, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	COA4 (1-87, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMSTSVPO GHTWTQRVKK DDEEEDPLDQ LISRSGCAAS HFAVQECMAQ HQDWRQCQPQ VQAFKDCMSE QQARRQEELQ RRQEQA GAHH
Tag:	His-tag
Predicted MW:	12.5 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human COA4 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_057649
Locus ID:	51287
UniProt ID:	Q9NYJ1 , A0A024R5J6
Cytogenetics:	11q13.4
Synonyms:	CHCHD8; CMC3; E2IG2
Summary:	Putative COX assembly factor.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

