

Product datasheet for **AR51549PU-S**

ABHD12B (1-255, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	ABHD12B (1-255, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMLGIWHT VPSCRGEDAK GKDCCWYEEA LRDGNPIIVY LHGSAEHRAA SHRLKLVKVL SDGGFHVLSV DYRGFGDSTG KPTEEGLTTD AICVYEWTKA RSGITPVCLW GHSLGTGVAT NAAKVLEEK CPVDAIVLEA PFTNMWVASI NYPLLLKIYRN IPGFLRTLMD ALRKDKIIFP NDENVKFLSS PLLILHGEDD RTVPLEYGKK LYEIARNAYR NKERVKMFVIF PPGFQHNLCC KSPTLLITVR DFLSKQWS
Tag:	His-tag
Predicted MW:	31.0 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human ABHD12B protein, fused to His-tag at N-terminus, was expressed in E.coli.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001193602
Locus ID:	145447
UniProt ID:	Q7Z5M8
Cytogenetics:	14q22.1
Synonyms:	BEM46L3; C14orf29; c14_5314



[View online »](#)

Summary: ABHD12B belongs to the serine esterase family.

Product images:

