

Product datasheet for **AR51562PU-S**

THAP11 (1-314, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	THAP11 (1-314, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMPGFTCC VPGCYNNSHR DKALHFYTFP KDAELRRLWL KNVSRAGVSG CFSTFQPTTG HRLCSVHFQG GRKTYTVRVP TIFPLRGVNE RKVARRPAGA AAARRRQQQQ QQQQQQQQQQ QQQQQQQQQQ QQQQQSSPSA STAQTAQLQP NLVSASAAVL LTLQATVDSS QAPGSVQPAP ITPTGEDVKP IDLTVQVEFA AAEGAAAAAA ASELQAATAG LEAAECPMGP QLVVGEEGF PDTGSDHSYS LSSGTTEEEL LRKLNEQRDI LALMEVKMKE MKGSIRHLRL TEAKLREELR EKDRLLAMAV IRKKHGM
Tag:	His-tag
Predicted MW:	36.8 kDa
Concentration:	lot specific
Purity:	>80% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human THAP11 protein, fused to His-tag at N-terminus, was expressed in E.coli.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_065190
Locus ID:	57215
UniProt ID:	Q96EK4
Cytogenetics:	16q22.1
Synonyms:	CTG-B43a; CTG-B45d; HRIHFB2206; RONIN



[View online »](#)

Summary:

The protein encoded by this gene contains a THAP domain, which is a conserved DNA-binding domain that has striking similarity to the site-specific DNA-binding domain (DBD) of *Drosophila* P element transposases. [provided by RefSeq, Jul 2008]

Product images: