

Product datasheet for AR51563PU-S

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200 Rockville, MD 20850 LIS

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

BMP8B (264-402, His-tag) Human Protein

Product data:

Product Type: Recombinant Proteins

Description: BMP8B (264-402, His-tag) human recombinant protein, 0.1 mg

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

MGSSHHHHHH SSGLVPRGSH MGSAVRPLRR RQPKKSNELP QANRLPGIFD DVHGSHGRQV CRRHELYVSF QDLGWLDWVI APQGYSAYYC EGECSFPLDS CMNATNHAIL QSLVHLMMPD

CRRHELYVSF QDLGWLDWVI APQGYSAYYC EGECSFPLDS CIVINATINHAIL QSLVHLIVIIVIPL

AVPKACCAPT KLSATSVLYY DSSNNVILRK HRNMVVKACG CH

Tag: His-tag
Predicted MW: 18.1 kDa
Concentration: lot specific

Purity: >90% by SDS - PAGE

Buffer: Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% glycerol

Preparation: Liquid purified protein

Protein Description: Recombinant human BMP8B protein, fused to His-tag at N-terminus, was expressed in E.coli.

Storage: Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid

repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: <u>NP 001711</u>

Locus ID: 656

UniProt ID: P34820
Cytogenetics: 1p34.2

Synonyms: BMP8; OP2





Summary: This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta)

superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to

recruitment and activation of SMAD family transcription factors that regulate gene

expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer. The encoded protein stimulates thermogenesis in brown adipose tissue. Expression of this gene may be downregulated in pancreatic cancer. This gene

may have arose from a gene duplication event and its gene duplicate is also present on

chromosome 1. [provided by RefSeq, Jul 2016]

Protein Families: Adult stem cells, Cancer stem cells, Druggable Genome, Embryonic stem cells, ES Cell

Differentiation/IPS, Induced pluripotent stem cells, Secreted Protein, Stem cell relevant

signaling - TGFb/BMP signaling pathway

Protein Pathways: Hedgehog signaling pathway, TGF-beta signaling pathway

Product images:

