

## Product datasheet for **AR51659PU-N**

### MRPL48 (29-212, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	MRPL48 (29-212, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSSGEKPIY SVGGILLSIS RPYKTKPTHG IGKYKHLIKA EEPKKKKGKV EVRAINLGTD YEYGVLNHL TAYDMTLAES YAQYVHNLN SLSIKVEESY AMPTKTIEVL QLQDQGSKML LDSVLTTHR VVQISGLSAT FAEIFLEIIQ SSLPEGVRLS VKEHTEEDFK GRFKARPELE ELLAKLK
Tag:	His-tag
Predicted MW:	23.1 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol.
Preparation:	Liquid purified protein
Protein Description:	Recombinant human MRPL48 protein, fused to His-tag at N-terminus, was expressed in E.coli.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_001305427</a>
Locus ID:	51642
UniProt ID:	<a href="#">Q96GC5</a>
Cytogenetics:	11q13.4
Synonyms:	CGI-118; HSPC290; L48MT; MRP-L48



[View online »](#)

**Summary:**

Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. A pseudogene corresponding to this gene is found on chromosome 6p. Several transcript variants, some protein-coding and some non-protein coding, have been found for this gene. [provided by RefSeq, Jan 2016]

**Product images:**