

Product datasheet for AR51691PU-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

N6AMT1 / HEMK2 (1-214, His-tag) Human Protein

Product data:

Product Type: Recombinant Proteins

Description: N6AMT1 / HEMK2 (1-214, His-tag) human recombinant protein, 0.1 mg

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

MGSSHHHHHH SSGLVPRGSH MGSMAGENFA TPFHGHVGRG AFSDVYEPAE DTFLLLNALE

AAAAELAGVE ICLEVGSGSG VVSAFLASMI GPQALYMCTD INPEAAACTL ETARCNKVHI QPVITDLVKG LLPRLTEKVD LLVFNPPYVV TPPQEVGSHG IEAAWAGGRN GREVMDRFFP

LVPDLLSPRG LFYLVTIKEN NPEEILKIMK TKGLQGTTAL SRQAGQETLS VLKFTKS

Tag: His-tag

Predicted MW: 25.3 kDa

Concentration: lot specific

Purity: >85% by SDS - PAGE

Buffer: Presentation State: Purified

Buffer System: Liquid, In Phosphate buffered saline (pH 7.4) containing 20% glycerol, 1 mM

DTT

Protein Description: Recombinant human N6ATM1 protein, fused to His-tag at N-terminus, was expressed in E.coli

and purified by using conventional chromatography techniques.

Storage: Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid

repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: NP 037372

 Locus ID:
 29104

 UniProt ID:
 Q9Y5N5

 Cytogenetics:
 21q21.3

Synonyms: C21orf127; HEMK2; KMT9; m.HsaHemK2P; MTQ2; N6AMT; PRED28; PrmC





Summary:

This gene encodes an N(6)-adenine-specific DNA methyltransferase. The encoded enzyme may be involved in the methylation of release factor I during translation termination. This enzyme is also involved in converting the arsenic metabolite monomethylarsonous acid to the less toxic dimethylarsonic acid. Alternative splicing pf this gene results in multiple transcript variants. A related pseudogene has been identified on chromosome 11. [provided by RefSeq, Jul 2014]

Protein Families:

Druggable Genome

Product images:

