

## Product datasheet for **AR51701PU-S**

### HDGFL1 (1-251, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	HDGFL1 (1-251, His-tag) human recombinant protein, 20 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMSAYGMP MYKSGDLVFA KLKGYAHWPA RIEHMTQPNR YQVFFFGTHE TAFLSPKRLF PYKECKEKFG KPNKRRGFSA GLWEIENNPT VQASDCPLAS EKSGSDGPWP EPEAAEGDED KPTHAGGGGD ELGKPDDDKP TEEKGPLKR SAGDPPEDAP KRPKEAAPDQ EEEAEAERAA EAERAAAAAA ATAVDEESPF LVAVENGSA SEPGLVCEPP QPEEEELREE EVADEEASQE WHAEAPGGGD RDSL
Tag:	His-tag
Predicted MW:	29.6 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified Buffer System: Liquid, In Phosphate buffered saline (pH 7.4) containing 10% glycerol
Protein Description:	Recombinant human HDGFL1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_612641</a>
Locus ID:	154150
UniProt ID:	<a href="#">Q5TGJ6</a> , <a href="#">A0A140VJK8</a>
Cytogenetics:	6p22.3
Synonyms:	dj309H15.1; Hdgfrp1; HRP-1; PWWP1
Summary:	HDGFL1 belongs to the HDGF family and contains 1 PWWP domain.



[View online »](#)

**Product images:**

