

Product datasheet for **AR51731PU-S**

BRD3 (1-416, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	BRD3 (1-416, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMSTATTV APAGIPATPG PVNPPPPEVS NPSKPGRKTN QLQYMQNVVV KTLWKHQFAW PFYQPVDAIK LNLDPYHKII KNPMMDMGTIK KRLENNYYWS ASECMQDFNT MFTNCYIYNK PTDDIVLMAQ ALEKIFLQKV AQMPQEEVEL LPPAPKGKGR KPAAGAQSAG TQQVAAVSSV SPATPFQSVPTVSQTPVIA ATPVPTITAN VTSVPVPPAA APPPPATPIV PVVPPTPPVW KKKGVKRRKAD TTTPTTSAIT ASRSESPPL SDPKQAKVVA RRESGGRPIK PPKKDLEDGE VPQHAGKKGK LSEHLRYCDS ILREMLSKKH AAYAWPFYKP VDAEALHLD YHDIHKHPMD LSTVKKRMDG REYPDAGGFA ADVRLMFSNC YKYNPPDHEV VAMARKLQDV FEMRFAKMP
Tag:	His-tag
Predicted MW:	48.1 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: PBS, pH 7.4 containing 10% Glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant Human BRD3 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_031397
Locus ID:	8019
UniProt ID:	Q15059 , A0A024R8H6
Cytogenetics:	9q34.2



[View online »](#)

Synonyms: ORFX; RING3L

Summary: This gene was identified based on its homology to the gene encoding the RING3 protein, a serine/threonine kinase. The gene localizes to 9q34, a region which contains several major histocompatibility complex (MHC) genes. The function of the encoded protein is not known. [provided by RefSeq, Jul 2008]

Protein Families: Protein Kinase

Product images:

