

Product datasheet for **AR51754PU-S**

TATDN1 (1-297, His-tag) Human Protein

Product data:


Product Type:	Recombinant Proteins
Description:	TATDN1 (1-297, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMSRFKFI DIGINLTDPM FRGIYRGVQK HQDDLQDVIG RAVEIGVKKF MITGGNLQDS KDALHLAQTN GMFFSTVGCH PTRCGEFEKN NPDLYLKELL NLAENNKGVKVAIGECGLDF DRLQFCPKDT QLKYFEKQFE LSEQTKLPMF LHCNRNSHAEF LDIMKRNDR CVGGVHVSFD GTKEAAAALI DLDLYIGFNG CSLKTEANLE VLKSIPSEKL MIETDAPWCG VKSTHAGSKY IRTAFPTKKK WESGHCLKDR NEPCHIIQIL EIMSAVRDED PLELANTLYN NTIKVFFPGI
Tag:	His-tag
Predicted MW:	36.0 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: Phosphate buffered saline (pH 7.4) containing, 10% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human TATDN1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001139632
Locus ID:	83940
UniProt ID:	Q6P1N9
Cytogenetics:	8q24.13
Synonyms:	CDA11



[View online »](#)

Summary: Putative deoxyribonuclease.[UniProtKB/Swiss-Prot Function]

Product images:



15% SDS-PAGE (3ug)