

Product datasheet for **AR51768PU-S**

RAB3IL1 (1-605, His-tag) Human Protein

Product data:


Product Type:	Recombinant Proteins
Description:	RAB3IL1 (1-605, His-tag) human recombinant protein, 20 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MMDLELPPPG LPSQQDMDLI DILWRQDIDL GVSREVFDFS QRRKEYELEK QKKLEKERQE QLQKEQEKAF FAQLQLDEET GEFLPIQPAQ HIQSETSGSA NYSQVAHIPK SDALYFDDCM QLLAQTFPFV DDNEVSSATF QSLVPDIPGH IESPVFIATN QAQSPETSVA QVAPVDLDGM QQDIEQWEE LLSIPELQCL NIENDKLVET TMVPSPEAKL TEVDNYHFYS SIPSMEKEVG NCSPHFLNAF EDSFSSILST EDPNQLTVNS LNSDATVNTD FGDEFYSAFI AEPSISNSMP SPATLSHLS ELLNGPIDVS DLSLCKAFNQ NHPESTAEFN DSDSGISLNT SPSVASPEHS VESSYGDTL LGLSDSEVEE LDSAPGSVKQ NGPKTPVHSS GDMVQPLSPS QGQSTHVHDA QCENTPEKEL PVSPGHRKTP FTKDKHSSRL EAHLTRDELAKALHIPFPV EKIINLPVVD FNEMMSKEQF NEAQLALIRD IRRRGKKNVA AQNCRKRKLE NIVELEQLD HLKDEKEKLL KEKGENDKSL HLLKKQLSTL YLEVFSMLRD EDGKPYSPSE YSLQQTRDGN VFLVPKSKKP DVKKN
Tag:	His-tag
Predicted MW:	69.9 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: Phosphate buffer saline (pH 7.4) containing 10% glycerol.
Preparation:	Liquid purified protein
Protein Description:	Recombinant human NRF2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001258615
Locus ID:	5866



[View online »](#)

UniProt ID:	Q8TBN0
Cytogenetics:	11q12.2-q12.3
Synonyms:	GRAB
Summary:	This gene encodes a guanine nucleotide exchange factor for the ras-related protein Rab3A. The encoded protein binds Rab3a and the inositol hexakisphosphate kinase InsP6K1. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 7. [provided by RefSeq, Nov 2012]

Product images:



15% SDS-PAGE (3ug)