

Product datasheet for **AR51775PU-S**

Interleukin-29 / IL29 (20-200, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Interleukin-29 / IL29 (20-200, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Tag:	His-tag
Predicted MW:	22.7 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: Liquid, In 20mM Tris-HCl (pH8.0) containing 10% glycerol.
Preparation:	Liquid purified protein
Protein Description:	Recombinant human IFNL1, fused to His-tag at N-terminus, was expressed in E.coli.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_742152
Locus ID:	282618
Cytogenetics:	19q13.2
Synonyms:	IL-29; IL29
Summary:	This gene encodes a cytokine distantly related to type I interferons and the IL-10 family. This gene, interleukin 28A (IL28A), and interleukin 28B (IL28B) are three closely related cytokine genes that form a cytokine gene cluster on a chromosomal region mapped to 19q13. Expression of the cytokines encoded by the three genes can be induced by viral infection. All three cytokines have been shown to interact with a heterodimeric class II cytokine receptor that consists of interleukin 10 receptor, beta (IL10RB) and interleukin 28 receptor, alpha (IL28RA). [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane



[View online »](#)

Protein Pathways: Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway

Product images:

