

Product datasheet for **AR51878PU-S**

Apolipoprotein A I / APO AI (25-264, His-tag) Mouse Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Apolipoprotein A I / APO AI (25-264, His-tag) mouse protein, 0.1 mg
Species:	Mouse
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSDEPQSQW DKVKDFANVY VDAVKDSGRD YVSQFESSLL GQQLNLLNLE NWDTLGSTVS QLQERLGPLT RDFWDNLEKE TDWVRQEMNK DLEEVKQKVQ PYLDEFQKKW KEDVELYRQK VAPLGAELQE SARQKLQELQ GRLSPVAEEF RDRMRTHVDS LRTQLAPHSE QMRESLAQRL AELKSNPTLN EYHTRAKTHL KTLGEEKARPA LEDLRHSLMP MLETLKTQVQ SVIDKASETL TAQ
Tag:	His-tag
Predicted MW:	30.3 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: Liquid, In Phosphate buffered saline (pH 7.4) containing 20% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_033822
Locus ID:	11806
UniProt ID:	Q00623
Cytogenetics:	9 25.36 cM
Synonyms:	AI; Alp-1; Ap; apo-AI; APOA-1; apoA-I; Brp-; Brp-14; Ltw-; Ltw-1; Lvtw; Lvtw-1; Se; Sep; Sep-1; Sep-2; Sep2



[View online »](#)

Summary:

This gene encodes a preproprotein that is proteolytically cleaved to yield a signal peptide and a proprotein that is subsequently processed to generate the active mature peptide. The encoded protein is the major protein component of plasma high density lipoprotein (HDL). This protein facilitates the removal of cholesterol and other fats from tissues by transporting them to the liver for excretion. This protein is a cofactor for lecithin cholesterolacyltransferase, an enzyme that catalyzes the conversion of free cholesterol to cholesteryl esters. Mutations in this gene in humans causes familial HDL deficiency, Tangier disease and familial visceral amyloidosis. Similar clinical features are exhibited by mice with mutations in this gene. This gene is clustered with three other apolipoprotein genes on chromosome 9. [provided by RefSeq, Dec 2013]

Product images: