

Product datasheet for **AR51953PU-S**

LRPAP1 / A2MRAP (29-360, His-tag) Mouse Protein

Product data:


Product Type:	Recombinant Proteins
Description:	LRPAP1 / A2MRAP (29-360, His-tag) mouse recombinant protein, 0.1 mg
Species:	Mouse
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSQPIAGHG GKYSREKNEP EMAAKRESGE EFRMEKLNQL WEKAKRLHLS PVRLAELHSD LKIQERDELN WKCLKVEGLD KDGEKEAKLI HNLNVILARY GLDGRKDAQM VHSNALNEDT QDELGDPRLE KLWHKAKTSG KFSSEELDKL WREFLHYKEK IQEYNVLLDT LSRAEEGYEN LLSPSDMAHI KSDTLISKHS ELKDRLRSIN QGLDRLRKVS HQGYGSTTEF EEPRLVDLWD LAQSANFTEK ELESFRELK HFEAKIEKHN HYQKQLEISH QCLKHVESIG DPEHISRNKE KYVLLLEETK ELGYKVKKHL QDLSSRVSRA RHNEL
Tag:	His-tag
Predicted MW:	41.4 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: Liquid, In Phosphate buffered saline (pH 7.4) containing 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant mouse Lspap1, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_038615
Locus ID:	16976
UniProt ID:	P55302 , Q52KI7
Cytogenetics:	5 18.01 cM
Synonyms:	AA617339; AI790446; AU042172; C77774; HBP44; RAP



[View online »](#)

Summary: Molecular chaperone for LDL receptor-related proteins that may regulate their ligand binding activity along the secretory pathway.[UniProtKB/Swiss-Prot Function]

Product images:



15% SDS-PAGE (3ug)