

Product datasheet for **AR51980PU-S**

Alpha-galactosidase A / GLA (32-429, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Alpha-galactosidase A / GLA (32-429, His-tag) human protein, 50 µg
Species:	Human
Expression Host:	Insect
Expression cDNA Clone or AA Sequence:	LDNGLARTPT MGWLHWERFM CNLDCQEEP SCISEKLFME MAELMVSEGW KDAGYEYLCI DDCWMAPQRD SEGRLQADPQ RFPHGIRQLA NYVHSKGLKL GIYADVGNKT CAGFPGSFGY YDIDAQTFAD WGVDLLKFDG CYCDSLENLA DGYKHMSLAL NRTGRSIVYS CEWPLYMWPF QKPNYTEIRQ YCNHWRNFAD IDDSWKSIS ILDWTSFNQE RIVDVAGPGG WNDPDMLVIG NFGLSWNQOV TQMALWAIMA APLFMSNDLR HISPQAKALL QDKDVAINQ DPLGKQGYQL RQGDNFEVWE RPLSGLAWAV AMINRQEIGG PRSYTIAVAS LGKGVACNPA CFITQLLPVK RKLGFYEWTS RLRSHINPTG TVLLQLENTM QMSLKDLLVE HHHHHH
Tag:	His-tag
Predicted MW:	46.4 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: Phosphate Buffered Saline (pH 7.4) containing 10% glycerol.
Endotoxin:	< 1.0 EU per 1 microgram of protein (determined by LAL method)
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_000160
Locus ID:	2717
UniProt ID:	P06280 , Q53Y83
Cytogenetics:	Xq22.1
Synonyms:	GALA



[View online »](#)

Summary:


This gene encodes a homodimeric glycoprotein that hydrolyses the terminal alpha-galactosyl moieties from glycolipids and glycoproteins. This enzyme predominantly hydrolyzes ceramide trihexoside, and it can catalyze the hydrolysis of melibiose into galactose and glucose. A variety of mutations in this gene affect the synthesis, processing, and stability of this enzyme, which causes Fabry disease, a rare lysosomal storage disorder that results from a failure to catabolize alpha-D-galactosyl glycolipid moieties. [provided by RefSeq, Jul 2008]

Protein Families:

Druggable Genome

Protein Pathways:

Galactose metabolism, Glycerolipid metabolism, Glycosphingolipid biosynthesis - globo series, Lysosome, Sphingolipid metabolism

Product images:

15% SDS-PAGE (3ug)