

## Product datasheet for **AR52001PU-S**

### CD106 / VCAM1 (25-698, His-tag) Mouse Protein

#### Product data:

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	CD106 / VCAM1 (25-698, His-tag) mouse protein, 50 µg
<b>Species:</b>	Mouse
<b>Expression Host:</b>	Insect
<b>Expression cDNA Clone or AA Sequence:</b>	FKIEISPEYK TIAQIGDSMA LTCSTTGCEs PLFSWRTQID SPLNAKVRTE GSKSVLTMEP VSFENEHSYL CTATCGSGKL ERSIHVDIYS FPKDPEIQFS GPLEVGGKPVV VKCLAPDIYP VYRLEIDLK GDQLMNRQEF SSEEMTKSLE TKSLEVTFTP VIEDIGKALV CRAKLHIDQI DSTLKERETV KELQVYISPR NTTISVHPST RLQEGGAVTM TCSSEGLPAP EIFWGRKLDN EVLQLLSGNA TLTIAMRME DSGVVVCEGV NLIGRDKAEV ELVQEKPFV VDISPGSQVA AQVGDSVLT CAAIGCDSPS FSWRTQTDSP LNGVVRNEGA KSTLVLSVSG FEDEHSYLCA VTCLQRTLEK RTQVEVYFSP EDPVIKMSGP LVHGRPVTVN CTVPNVYFPD HLEIELLKGE TTLMKKYFLE EMGIKSLETK ILETTFIPTI EDTGKSLVCL ARLHSGEMES EPKQRQSVQP LYVNVAPKET TIWVSPSPIL EEGSPVNLTC SSDGIPAPKI LWSRQLNNGE LQPLSENTTL TFMSTKRDDS GIYVCEGINE AGISRKSVEL IIQVSPKDIQ LTVFPSKSVK EGDTVIISCT CGNVPETWII LKKKAKTGDM VLKSVGDGSYT IRQAQLQDAG IYECESKTEV GSQLRSLTLD VKGKEHNKNY FSPELEHHHH HH
<b>Tag:</b>	His-tag
<b>Predicted MW:</b>	75.4 kDa
<b>Concentration:</b>	lot specific
<b>Purity:</b>	>95% by SDS – PAGE
<b>Buffer:</b>	Presentation State: Purified State: Liquid purified protein Buffer System: Phosphate Buffered Saline (pH 7.4) containing 10% glycerol.
<b>Endotoxin:</b>	< 1.0 EU per 1 microgram of protein (determined by LAL method)
<b>Preparation:</b>	Liquid purified protein
<b>Storage:</b>	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.
<b>RefSeq:</b>	<a href="#">NP_035823</a>
<b>Locus ID:</b>	22329



[View online »](#)

UniProt ID:	<a href="#">Q3UPN1</a>
Cytogenetics:	3 50.17 cM
Synonyms:	CD106; Vcam-1
Summary:	VCAM1, also known as vascular cell adhesion protein 1, is a receptor for encephalomyocarditis virus on murine vascular endothelial cells. It has been identified as a receptor for the D variant of encephalomyocarditis (EMC-D) virus on vascular endothelial cells from the heart. Recombinant mouse VCAM1, fused to His-tag at C-terminus, was expressed in insect cell and purified by using conventional chromatography techniques.

### Product images:

