

Product datasheet for **AR52008PU-N**

CD172a / SIRPA (32-373, His-tag) Rat Protein

Product data:


Product Type:	Recombinant Proteins
Description:	CD172a / SIRPA (32-373, His-tag) rat protein, 0.25 mg
Species:	Rat
Expression Host:	Insect
Expression cDNA Clone or AA Sequence:	KELKVTQADK SVSVAAGDSA TLNCTVSSLT PVGPIKWFKG EGQNRSPIYS FIGGEHFPRI TNVSDATKRN NMDFSICISN VTPEDAGTY CVKFQKGIVE PDTEIKSGGG TTTYVLAKPS SPEVSGPDSR GSPGQTVNFT CKSYGFSPRN ITLKWLKNGK ELSHLETTIS SKSNVSYNIS STVSVKLSPE DIHSRVICEV AHVTLEGRPL NGTANFSNII RVSPTLKITQ QPLTPASQVN LTCQVQKFYP KALQLNWLEN GNLSRTDKPE HFTDNRDGTY NYTSLFLVNS SAHREDVVFT CQVEHDSQPA ITENHTVRAF AHSSSGGSME TIPDNNAYYN WNVEHHHHHH
Tag:	His-tag
Predicted MW:	38.5 kDa
Concentration:	lot specific
Purity:	>95% by SDS – PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: Phosphate Buffered Saline (pH 7.4) containing 10% glycerol.
Endotoxin:	< 1.0 EU per 1 microgram of protein (determined by LAL method)
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_037148
Locus ID:	25528
UniProt ID:	P97710
Cytogenetics:	3q36
Synonyms:	Bit; Ptpns1; SHPS-1



[View online »](#)

Summary: plays a role in macrophage fusion [RGD, Feb 2006]

Product images:



15% SDS-PAGE (3ug)