

Product datasheet for **AR52059PU-S**

TPST2 (26-377, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	TPST2 (26-377, His-tag) human protein, 50 µg
Species:	Human
Expression cDNA Clone or AA Sequence:	ADPQQVLECR AVLAGLRSPR GAMRPEQEEL VMVGTNHVEY RYGKAMPLIF VGGVPRSGTT LMRAMLDAHP EVRCGEETRI IPRVLAMRQA WSKSGREKLR LDEAGVTDEV LDAAMQAFIL EVIAKHGEP A RVLCNKDPFT LKSSVYLSRL FPNSKFLLMV RDGRASVHSM ITRKVTIAGF DLSSYRDCLT KWNKAIEVMY AQCMEVGKEK CLPVYQEQLV LHPRRSLKLI LDFLGIAWSD AVLHHEDLIG KPGGVSLSKI ERSTDQVIK VNLEALSKWT GHIPGDVWRD MAQIAPMLAQ LGYDPYANPP NYGNPDPFVI NNTQRVLKGD YKTPANLKG Y FQVNQNSTSS HLGSSHHHHH H
Tag:	His-tag
Predicted MW:	40.4 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: Phosphate Buffered Saline (pH 7.4) containing 10% glycerol.
Endotoxin:	< 1.0 EU per 1 microgram of protein (determined by LAL method)
Preparation:	Liquid purified protein
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001008566
Locus ID:	8459
UniProt ID:	O60704 , A0A024R1G9
Cytogenetics:	22q12.1
Synonyms:	TANGO13B; TPST-2




[View online »](#)

Summary: The protein encoded by this gene catalyzes the O-sulfation of tyrosine residues within acidic regions of proteins. The encoded protein is a type II integral membrane protein found in the Golgi body. Alternative splicing produces multiple transcript variants encoding distinct isoforms. [provided by RefSeq, May 2018]

Protein Families: Transmembrane

Product images:



15% SDS-PAGE (3ug)