

Product datasheet for BA1000

RNA5-8SN2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Human Myeloma IgA1, Kappa. Single arc by IEP against antisera to whole Human serum, IgA1 and kappa. No reaction by IEP against IgA2, IgG, IgM, IgD, IgE and lambda.
Species:	Human
Expression Host:	Myeloma
Concentration:	2.81 mg/ml (OD280 nm, E1% = 13.4)
Purity:	> 95 % pure by SDS PAGE
Buffer:	10mM Sodium Phosphate, pH 7.4 containing 150mM Sodium Chloride
Preparation:	Liquid purified fraction



[View online »](#)