

Product datasheet for **BA1009S**

Desmin Chicken Protein

Product data:

Product Type:	Native Proteins
Description:	Desmin chicken protein, 0.1 mg
Species:	Chicken
Protein Source:	Stomach
Predicted MW:	53 kDa
Concentration:	lot specific
Purity:	>98% pure as determined by SDS-PAGE
Buffer:	Presentation State: Purified State: Lyophilized Buffer System: 10 mM Sodium Phosphate, pH 7.5, 6M urea, 2 mM DDT, 1 mM EDTA
Reconstitution Method:	Restore with distilled water. BA1009S: 100 µl (final volume 100 µl). BA1009 : 250 µl (final volume 250 µl).
Preparation:	Lyophilized
Applications:	Protein standard in 1D and 2D SDS gel electrophoresis. Immunoassay. Antigen for immunization.
Protein Description:	Purified Chicken Desmin.
Note:	Isoelectric point: pHI 5.4
Storage:	Store lyophilized product at 2-8°C . After reconstitution store at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	XP_015145578
Locus ID:	395906
Summary:	desmin



[View online »](#)