

Product datasheet for BA1016

Cytokeratin 19 Human Protein

Product data:

Concentration:

Product Type: Recombinant Proteins

Description: Cytokeratin 19 human recombinant protein, 0.25 mg

Species:HumanExpression Host:E. coliPredicted MW:44 kDa

Purity: >95% determined by SDS gelelectrophoresis

lot specific

Buffer: Presentation State: Purified

State: Lyophilized purified protein

Reconstitution Method: BA1016: Restore with 175 μl distilled water (final volume 250 μl).

BA1016S: Restore with 70 μ l distilled water (final volume 100 μ l).

Final Solution: 30mM Tris/HCl pH 8, 9.5M Urea, 2mM DTT, 2mM EDTA, 10mM

Methylammonium Chloride.

Preparation: Lyophilized purified protein

Applications: Protein standard in 1D and 2D SDS gelelectrophoresis.

Immunoassays. Immunization.

Protocol: **Reconstitution to filaments** is performed by mixing equimolar amounts of keratins of type I and type II at concentrations of approx. 0.5 mg/ml, both dissolved in 9.5 M urea buffer (see above). Protofilaments and filament complexes are obtained by dialyzing the resulting polypeptide solution stepwise to a concentration of 4 M urea and then to low salt

condition (50 mM NaCl, 2 mM dithiothreitol, 10 mM Tris-HCl, pH 7.4).

For immunization purposes, the solution can be further dialyzed against PBS (phosphate

buffered saline, e.g. Dulbeccos PBS).

Protein Description: Recombinant Human keratin K19 (formerly also designated Cytokeratin 19).

Note: <u>Isoelectric Point</u>: pl 5.1



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Cytokeratin 19 Human Protein - BA1016

Storage: Store lyophilized at 2-8°C for 6 months or at -20°C long term.

After reconstitution store undiluted at 2-8°C for one month

or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: <u>NP 002267</u>

 Locus ID:
 3880

 UniProt ID:
 P08727

 Cytogenetics:
 17q21.2

Synonyms: CK19; K1CS; K19

Summary: The protein encoded by this gene is a member of the keratin family. The keratins are

intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. The type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. Unlike its related family members, this smallest known acidic cytokeratin is not paired with a basic cytokeratin in epithelial cells. It is specifically expressed in the periderm, the transiently superficial layer that envelopes the developing epidermis. The type I cytokeratins are clustered in a region of

chromosome 17q12-q21. [provided by RefSeq, Jul 2008]